

See discussions, stats, and author profiles for this publication at: <https://www.researchgate.net/publication/47297501>

Optimized Workflow for Preparation of APTS-Labeled N-Glycans Allowing High-Throughput Analysis of Human Plasma Glycomes using 48-Channel Multiplexed CGE-LIF

ARTICLE in JOURNAL OF PROTEOME RESEARCH · OCTOBER 2010

Impact Factor: 4.25 · DOI: 10.1021/pr100802f · Source: PubMed

CITATIONS

49

READS

75

8 AUTHORS, INCLUDING:



Lucia Renee Ruhaak

Leiden University Medical Centre

48 PUBLICATIONS 1,162 CITATIONS

SEE PROFILE



Erdmann Rapp

Max Planck Institute for Dynamics of Comp...

85 PUBLICATIONS 1,180 CITATIONS

SEE PROFILE



Manfred Wuhrer

Leiden University Medical Centre

195 PUBLICATIONS 5,320 CITATIONS

SEE PROFILE

Optimized Workflow for Preparation of APTS-Labeled N-Glycans Allowing High-Throughput Analysis of Human Plasma Glycomes using 48-Channel Multiplexed CGE-LIF

L. Renee Ruhaak,^{*,†} René Hennig,^{†,‡} Carolin Huhn,^{†,§} Matthias Borowiak,^{†,‡} Radboud J. E. M. Dolhain,^{||} André M. Deelder,[†] Erdmann Rapp,[‡] and Manfred Wuhrer[†]

Department of Parasitology, Biomolecular Mass Spectrometry Unit, Leiden University Medical Center, Leiden, The Netherlands, Max Planck Institute for Dynamics of Complex Technical Systems, Magdeburg, Germany, Central Division of Analytical Chemistry, Forschungszentrum Jülich, Jülich, Germany, and Department of Rheumatology, Erasmus MC University Medical Center, Rotterdam, The Netherlands

Received August 3, 2010

High-throughput methods for oligosaccharide analysis are required when searching for glycan-based biomarkers. Next to mass spectrometry-based methods, which allow fast and reproducible analysis of such compounds, further separation-based techniques are needed, which allow for quantitative analysis. Here, an optimized sample preparation method for N-glycan-profiling by multiplexed capillary gel electrophoresis with laser-induced fluorescence detection (CGE-LIF) was developed, enabling high-throughput glycosylation analysis. First, glycans are released enzymatically from denatured plasma glycoproteins. Second, glycans are labeled with APTS using 2-picoline borane as a nontoxic and efficient reducing agent. Reaction conditions are optimized for a high labeling efficiency, short handling times, and only limited loss of sialic acids. Third, samples are subjected to hydrophilic interaction chromatography (HILIC) purification at the 96-well plate format. Subsequently, purified APTS-labeled N-glycans are analyzed by CGE-LIF using a 48-capillary DNA sequencer. The method was found to be robust and suitable for high-throughput glycan analysis. Even though the method comprises two overnight incubations, 96 samples can be analyzed with an overall labor allocation time of 2.5 h. The method was applied to serum samples from a pregnant woman, which were sampled during first, second, and third trimesters of pregnancy, as well as 6 weeks, 3 months, and 6 months postpartum. Alterations in the glycosylation patterns were observed with gestation and time after delivery.

Keywords: glycan profiling • high-throughput N-glycan analysis • CGE-LIF • APTS • HILIC-SPE • plasma

Introduction

A large proportion of the human proteome is N-glycosylated.¹ N-glycans play a role in protein folding as well as protein solubility and are often crucial for the activity and function of the protein.² Therefore, analysis of protein glycosylation is gaining interest from biotechnological and pharmaceutical industry as well as clinical research. Applications in both fields require high-throughput analytical methods that allow fast and robust profiling of protein glycosylation on large sample sets.

Several strategies for fast and high-throughput analysis of N-glycans have been reported during the last years (e.g., refs 3–7). Most of them are based on liquid chromatography (LC) or capillary electrophoresis (CE) separations, which are typically performed on one sample at a time, making the analysis of

multiple samples time-consuming. Run times for fast LC glycan analyses are in the range of 30 min.^{4,5} N-glycan separations by chip-CE with run times of 3 min per sample were already demonstrated;⁷ however, to our knowledge, the hardware is not yet commercially available. A faster approach is direct mass spectrometry using MALDI ionization, where many samples can be spotted at once for subsequent sequential analysis, with analysis times of approximately 20 s per sample. However, no separation of isobaric structures can be obtained and, more importantly, quantitation using mass spectrometry is rather demanding in terms of expertise and equipment.^{6,8}

Multiplexed capillary gel electrophoresis with laser-induced fluorescence detection (CGE-LIF) using a DNA sequencer has been reported for the rapid analysis of 8-aminopyrene-1,3,6-trisulfonic acid (APTS)-labeled glycans.^{9–13} CGE-LIF experiments can be performed on up to 96 samples in parallel, depending on the type of instrument. Sample cleanup after APTS labeling is either omitted¹³ or being performed by size-exclusion chromatography. As the established protocols are rather tedious¹¹ or inefficient in removing the excess label,¹⁰ we were looking for efficient alternative approaches.

* To whom correspondence should be addressed. L. R. Ruhaak, Leiden University Medical Center, Department of Parasitology, Postal Zone L4-Q Albinusdreef 2 2333 ZA Leiden, The Netherlands. Tel: +31-71-5265077. Fax: +31-71-5266907. E-mail: lrruhaak@gmail.com.

[†] Leiden University Medical Center.

[‡] Max Planck Institute for Dynamics of Complex Technical Systems.

[§] Forschungszentrum Jülich.

^{||} Erasmus MC University Medical Center.

We here describe the development of a robust and facile high-throughput sample preparation method at the 96 well plate format for analysis of the total plasma N-glycome by CGE-LIF of APTS-labeled glycans. The efficacy of the labeling of plasma derived N-glycans with APTS was optimized using the recently reported nontoxic reducing agent 2-picoline borane.¹⁴ A polyacrylamide-based stationary phase (Biogel P10) was used in hydrophilic interaction chromatography (HILIC) mode to purify APTS-labeled glycans. Dimethylsulfoxide has been selected as an optimal solvent for repeatable and efficient CGE electrokinetic sample injection. Peak annotation of the total plasma N-glycome was achieved by HPLC fractionation and mass spectrometric analysis of the glycans prior to APTS-labeling and CGE-LIF analysis.

Materials and Methods

Materials. Dimethylsulfoxide (DMSO), Nonidet P-40 (NP-40), triethylamine (TEA), aminopyrene-1,3,6-trisulfonic acid (APTS), sodium cyanoborohydride (NaBH_3CN), 2-picoline borane, trifluoroacetic acid (TFA), and ammonium formate were obtained from Sigma-Aldrich (Zwijndrecht, The Netherlands). Sodium dodecyl sulfate (SDS) was bought from United States Biochemicals (Cleveland, OH). Glycan standards Glycan 18 (pentamannosidic N-glycan; MC0731), Glycan 27 (asialo, agalacto, bisected, core-fucosylated biantennary N-glycan; C0980), and Glycans 28 and 30 (asialo, monogalactosylated biantennary N-glycan; C0870) were purchased from Dextra Laboratories (Reading, U.K.). Glycan 17 (monosialylated, core-fucosylated biantennary N-glycan; M3800) was from Sigma/Aldrich. Glycan 34 (asialo, core-fucosylated biantennary N-glycan; GTP 0N-2A*F) was from TheraProteins (Oeiras, Portugal). PNGase F was obtained from Roche Diagnostics (Mannheim, Germany). Biogel P-10 was obtained from Bio-Rad (Veenendaal, The Netherlands), while glacial acetic acid, citric acid, microcrystalline cellulose and ethanol were from Merck (Darmstadt, Germany). Acetonitrile was purchased from Biosolve (Valkenswaard, The Netherlands). Sepharose CL-4B beads (45–165 μm) were from GE Healthcare (Uppsala, Sweden). Diol-bonded silica (Chromabond, Sorbenz, OH) was obtained from Macherey-Nagel, Düren, Germany. Hi-Di Formamide was from Applied Biosystems/Life Technologies (Carlsbad, CA). MQ (Milli-Q deionized water $R > 18.2 \text{ M}\Omega \text{ cm}^{-1}$) was used throughout (Millipore, Amsterdam, The Netherlands). Human plasma for method development was obtained from a healthy donor. Serum from a healthy woman with an uncomplicated pregnancy was collected in the PARA study¹⁵ and was used for glycosylation analysis.

Preparation of Oligosaccharides. N-glycans from human plasma and serum were prepared as described previously.¹⁶ Shortly, proteins from 10 μL of plasma or serum were denatured after addition of 20 μL 2% SDS by incubation at 60 °C for 10 min. Subsequently, 10 μL 4% NP-40 and 0.5 mU of PNGase F in 10 μL 5 \times PBS were added to the samples. The samples were incubated overnight at 37 °C for N-glycan release.

Labeling of Oligosaccharides. Labeling of oligosaccharides was performed as published¹⁷ with slight modifications: 2 μL of N-glycan solution were mixed with 2 μL of a freshly prepared solution of label (APTS; 20 mM in 3.6 M citric acid) in a V-bottom polypropylene 96-well plate (Westburg, Leusden, The Netherlands). Two μL aliquots of freshly prepared reducing agent solution (different molarities of NaBH_3CN or 2-picoline borane¹⁴ in DMSO) were added, the plate was sealed using adhesive tape and after 5 min of shaking, the samples were

incubated at 37 °C for 16 h. To stop the reaction, 50 μL of acetonitrile/water (80:20 v/v) were added and the samples were mixed for 5 min.

HILIC-Solid Phase Extraction (SPE). Free label and reducing agent were removed from the samples using HILIC-SPE. An amount of 100 μL of a 100 mg/mL suspension of cellulose, Sepharose, diol-bonded silica beads or Biogel P-10 in water/ethanol/acetonitrile (70:20:10, v/v) was applied to each well of a 0.45 μm GHP filter plate (Pall Corporation, Ann Arbor, MI). Solvent was removed by application of vacuum using a vacuum manifold (Millipore, Billerica, MA). All wells were prewashed using 5 \times 200 μL water, followed by equilibration using 3 \times 200 μL acetonitrile/water (80:20, v/v). The samples were loaded to the wells, and the plate was shaken for 5 min on a shaker to enhance glycan binding. The wells were subsequently washed using 5 \times 200 μL acetonitrile/100 mM triethylamine (TEA) adjusted to pH 8.5 with acetic acid (80:20, v/v), followed by 3 \times 200 μL acetonitrile/water (80:20, v/v). Washing steps were performed by addition of solutions, incubation for 30 s, and removal of solvent by vacuum. The incubation steps were observed to result in better yields and repeatability. Water (100 μL) was applied followed by a 5 min incubation on the shaker (to allow swelling of the Biogel P-10 particles). Thereafter, 200 μL water were added followed by a 5 min incubation on the shaker and collection of eluates by vacuum in a 96-well V-bottom polypropylene deep well plate (Westburg, Leusden, The Netherlands). Another 200 μL water were added, followed by a 5 min incubation on the shaker and elution into the deep-well plate. The combined eluates were either analyzed immediately by CGE-LIF or stored at –20 °C until usage.

CGE-LIF using ABI-3730 DNA Sequencing Equipment. Two μL of N-glycan eluate were added to 60 μL of water, DMSO or Hi-Di Formamide in a PCR plate (Thermo Fischer Scientific via Westburg, Leusden, The Netherlands). Plates were sealed and turned several times for thorough mixing and subsequently centrifuged prior to analysis using an ABI-3730 DNA sequencer (Applied Biosystems) to avoid air bubbles at the bottom of the wells. The injection voltage was set to 7.5 kV, while the running voltage was 10 kV. The system was equipped with a 48 channel array with capillaries of 50 cm in length, and the capillaries were filled with POP-7 buffer (Applied Biosystems). The 3730 running buffer was obtained from Applied Biosystems. For peak annotation GeneScan 500 ROX size standard (Applied Biosystems, 2 μL of a 1:100 dilution) was added as the internal standard prior to analysis of each sample. Data was collected with a frequency of 10 Hz for 50 min.

Data Processing of DNA Sequencer Data. Data files were converted to .xml files using DataFileConverter, which is supplied by Applied Biosystems and then loaded into the in-house Matlab-based data processing tool “glyXtool” (developed at the Max Planck Institute for Dynamics of Complex Technical Systems, Magdeburg on Matlab v7.1, The Mathworks, Inc., Natick, MA). After smoothing and background adjustment, peak heights were normalized to the sum of all peak heights in the glycan profile.

Porous Graphitic Carbon (PGC)- SPE of Unlabeled Glycans. All 96 N-glycan releases from a 96 well plate were pooled. Each sample of the plate represented an N-glycan release from 10 μL of a human plasma standard sample. The pooled N-glycans were purified by PGC according to.¹⁸ Briefly, a PGC cartridge (Alltech, Ridderkerk, The Netherlands) was first washed using 2 \times 3 mL of acetonitrile/water (80:20, v/v), then conditioned using 3 \times 3 mL of water. Sample was applied in

High-Throughput APTS-Labeled N-Glycan Preparation

water and the cartridge was washed using 4×3 mL 0.1% TFA in water. N-glycans were eluted using 3×1 mL of acetonitrile/water (50:50, v/v) containing 0.1% TFA. The eluate was brought to dryness and reconstituted in 1 mL acetonitrile/water (80:20, v/v) prior to fractionation via HILIC-HPLC.

HILIC-HPLC Separation of Unlabeled Glycans. N-glycans were fractionated by HILIC-HPLC (TSK amide 80, $3 \mu\text{m}$, 150 mm \times 4.6 mm inner diameter column; Tosoh Bioscience, Stuttgart, Germany). The column was run at 300 $\mu\text{L}/\text{min}$, and a binary gradient was applied using acetonitrile as solvent A, and 50 mM ammonium formate pH 4.4 as solvent B. The following gradient conditions were applied: 0 min 30% solvent B; 5 min 30% solvent B; 50 min 50% solvent B; 51 min 80% solvent B; 54 min 80% solvent B; 55 min 30% solvent B. Fraction collection started at 35 min, and 20 s fractions were collected in a 96-deep well plate. Fractions were dried by vacuum centrifugation.

MALDI-FTICR-MS. HILIC-HPLC fractions were reconstituted in 200 μL of water, and 1 μL of each fraction was spotted on an anchorchip MALDI target plate (Bruker Daltonics, Bremen, Germany). Spots were overlaid with 1 μL DHB (from Bruker; 20 mg/mL in 50% acetonitrile) and left to dry. MALDI-FTICR-MS experiments were performed in the positive-ion mode on a 9.4 T FTICR APEX-ultra mass spectrometer equipped with a dual ESI/MALDI ion source (Apollo II), a quadrupole mass filter and a smartbeam laser system (Bruker Daltonics). The intermediate pressure MALDI source enables the detection of glycoconjugates containing labile monosaccharides.¹⁹ All experiments used a laser spot size of approximately 150 μm and a laser repetition rate of 200 Hz. The quadrupole was operated in RF-only mode.

A customized pulse program was used to perform the MALDI-FTICR-MS experiments, as previously described.¹⁹ Briefly, ions from 50 laser shots were accumulated in a hexapole and then transferred to the collision cell. This cycle was performed nine times, accumulating ions from 450 laser shots in the collision cell. The accumulated ions were then transferred to the ICR cell for a mass analysis scan. Each spectrum is the sum of eight scans. Data acquisition was performed using ApexControl 3.0.0 expert software (Bruker Daltonics), while DataAnalysis V4.0 (Bruker Daltonics) was used for data processing.

Results

Total Plasma N-Glycome Analysis by CGE-LIF. Total plasma N-glycome profiles were analyzed by CGE-LIF using a 48 channel DNA sequencer after N-glycan release by PNGase F, APTS-labeling, and sample purification by gel filtration or HILIC-SPE (Figure 1). We here describe the optimization of the sample preparation and analysis. Mild labeling conditions were chosen which to a large extent preserve the sialic acid residues. A polyacrylamide-based stationary phase allowed the fast purification of the APTS-labeled glycans at the 96 well plate format. Injection conditions were optimized to maximize injection of APTS-labeled glycans. Finally, preparative HILIC-HPLC and MALDI-FTICR-MS were applied for peak assignment. A typical electropherogram showing the obtained total plasma N-glycome profile is shown in Figure 2. In the migration time range between 18 and 32 min, more than 30 peaks were observed, most of them baseline separated due to the large resolving power of CGE. To judge the optimization steps, we analyzed the electropherograms on three different levels:




Figure 1. Schematic overview of the optimized sample preparation procedure for APTS-labeled N-glycans from plasma.




Figure 2. Typical electropherogram of plasma derived N-glycans labeled with APTS according to the optimized procedure. Peaks are annotated corresponding to Table 1.

(1) The signal intensity of Peak 2 was monitored as an indicator of the overall glycan yield.

(2) To preclude a possible bias of the sample preparation method for sialylated/nonsialylated glycans and to register the possible loss of sialic acids, the ratio of Peak 2 (disialylated) to Peak 32 (nonsialylated counterpart) was determined (ratio Peak 2/Peak 32), see Table 1.

(3) We monitored the intensity of the APTS-peak in the electropherograms to judge the efficacy of sample purification in removing excess APTS. As this peak migrates fast, it appears in the electropherogram outside the N-glycan migration range shown in Figure 2.

Identification of Most Abundant Glycans. To achieve an assignment of the peaks observed in the total plasma N-glycan profile (Figure 2), their elution positions were compared with those of APTS-labeled standard glycans of known structure (Table 1). This was performed after a correction for shifts in migration times: the time scale was normalized using a DNA-based (basepair) size standard, which fluoresces at a different wavelength than APTS and therefore does not interfere with the glycan signals. Seven of the observed plasma glycan peaks could be annotated using commercially available standards. Moreover, a second, mass-spectrometry based strategy was used to compare the glycan signals observed in our study with the well-described composition of the human plasma N-glycome.^{4,10,11,13,16,20,21} To this end, N-glycans were released enzymatically from plasma glycoproteins, and the unlabeled glycans were preparatively separated using HILIC-HPLC. Ninety-

Table 1. Identification and Annotation of 35 N-Glycan Peaks from Human Plasma Proteins^a

Peak No.	Peak position in base pairs	Identified by spiking in a standard oligosaccharide	Measured N-glycan mass [M + Na] ⁺	Theoretical N-glycan mass [M + Na] ⁺	N-Glycan composition	N-Glycan structure	Present in LC-fractions analyzed by MS	Present in LC-fractions analyzed by CGE-LIF	Peak No.	Peak position in base pairs	Identified by spiking in a standard oligosaccharide	Measured N-glycan mass [M + Na] ⁺	Theoretical N-glycan mass [M + Na] ⁺	N-Glycan composition	N-Glycan structure	Present in LC-fractions analyzed by MS	Present in LC-fractions analyzed by CGE-LIF
1	172.2 (+/- 0.2)		2901.984	2902.000	H ₆ N ₃ S ₃		66 - 71	66 - 72	19	255.3 (+/- 0.3)		3122.097	3122.095	H ₇ N ₆ S ₂ F ₁		75 - 78	73 - 78
2	175.8 (+/- 0.3)		2245.765	2245.772	H ₄ N ₄ S ₂		56 - 66	56 - 67	20	256.2 (+/- 0.3)		2303.815	2303.814	H ₇ N ₅ S ₁ F ₁		54 - 56, 61 - 62	52 - 57, 60 - 63
3	177.9 (+/- 0.2)		3048.081	3048.058	H ₆ N ₅ S ₃ F ₁		69 - 76	69 - 77	21	257.4 (+/- 0.4)		2100.729	2100.735	H ₇ N ₄ S ₂ F ₁		48 - 60	48 - 52
4	181.6 (+/- 0.3)		2245.762	2245.772	H ₄ N ₄ S ₂		53 - 56	52 - 56	22	277.4 (+/- 0.4)		1419.483	1419.476	H ₄ N ₂		44 - 48	45 - 49
5	187.6 (+/- 0.3)		2391.822	2391.830	H ₅ N ₄ S ₂ F ₁		58 - 63	58 - 63	23	278.3 (+/- 0.4)		2319.800	2319.809	H ₆ N ₅ S ₁		62 - 72	61 - 71
6	191.8 (+/- 0.2)		2594.907	2594.910	H ₇ N ₅ S ₂ F ₁		60 - 63	60 - 64	24	281.9 (+/- 0.4)		1485.544	1485.534	H ₇ N ₄ F ₁		25 - 29	25 - 36
7	201.0 (+/- 0.3)		3267.167	3267.132	H ₇ N ₆ S ₃		71, 73 - 76	70 - 78	25	292.8 (+/- 0.4)		2465.845	2465.867	H ₆ N ₅ S ₂ F ₁		71 - 77	71 - 77
8	206.6 (+/- 0.2)		3413.160	3413.190	H ₇ N ₆ S ₂ F ₁		75 - 76	74 - 78	26	297.0 (+/- 0.3)		2465.845	2465.867	H ₆ N ₅ S ₂ F ₁		71 - 77	71 - 77
9	212.1 (+/- 0.3)		2610.892	2610.904	H ₆ N ₅ S ₂		62 - 78	65 - 77	27	297.2 (+/- 0.3)	X	1688.615	1688.613	H ₇ N ₃ F ₁		29 - 33	29 - 34
10	217.0 (+/- 0.3)		2610.892	2610.9044	H ₆ N ₅ S ₂		62 - 78	62 - 71	28	314.4 (+/- 0.3)	X	1647.593	1647.587	H ₆ N ₄ F ₁		36 - 41, 44, 46 - 48	34 - 45
11	223.6 (+/- 0.3)		2756.956	2756.962	H ₆ N ₅ S ₂ F ₁		67 - 76	66 - 77	29	318.5 (+/- 0.3)		2684.956	2684.941	H ₇ N ₆ S ₁		73 - 76	73 - 76
12	225.9 (+/- 0.4)		1954.669	1954.677	H ₇ N ₄ S ₁		44 - 65	44 - 49	30	320.8 (+/- 0.4)	X	1647.593	1647.587	H ₆ N ₄ F ₁		36 - 41, 44, 46 - 48	36 - 48
13	228.6 (+/- 0.4)		2756.956	2756.962	H ₆ N ₅ S ₂ F ₁		67 - 75	66 - 77	31	325.4 (+/- 0.2)		1850.645	1850.666	H ₆ N ₅ F ₁		38, 40	37 - 42
14	233.4 (+/- 0.3)		1954.669	1954.677	H ₇ N ₄ S ₁		44 - 65	49 - 70	32	327.2 (+/- 0.4)	X	1663.583	1663.581	H ₇ N ₄		44 - 46, 50 - 60	43 - 46, 50 - 60
15	239.7 (+/- 0.2)		1954.669	1954.677	H ₇ N ₄ S ₁		44 - 65	53 - 56	33	351.2 (+/- 0.2)		1905.630	1905.634	H ₆ N ₂		69 - 73	69 - 72
16	248.5 (+/- 0.2)		2976.049	2976.037	H ₇ N ₆ S ₂		71, 73, 75 - 76	70 - 76, 78	34	352.2 (+/- 0.4)	X	1809.644	1809.639	H ₇ N ₄ F ₁		45 - 56, 58 - 62	45 - 56
17	249.5 (+/- 0.3)	X	2100.729	2100.735	H ₇ N ₄ S ₁ F ₁		48 - 60	52 - 61, 63	35	361.8 (+/- 0.3)		2012.714	2012.719	H ₇ N ₅ F ₁		47 - 50, 54 - 56	47 - 50, 54 - 56
18	250.7 (+/- 0.3)	X	1257.444	1257.423	H ₅ N ₂		35 - 37	34 - 39									

^a Compositions are given in terms of hexose (H), N-acetylhexosamine (N), fucose (F) and N-acetyl neuraminic acid (S) and structural schemes are given in terms of galactose (yellow ball), N-acetyl glucosamine (blue square), mannose (green circle), fucose (red triangle) and sialic acid (purple diamond).

six fractions were collected and one aliquot of each fraction was analyzed by MALDI-FTICR-MS in the positive mode to obtain accurate masses. A second aliquot of each fraction was labeled with APTS and analyzed using CGE-LIF according to the optimized procedure. Mass spectrometric data of the unlabeled glycans as well as elution positions in HILIC-HPLC¹⁶ and previously published annotations of CGE-LIF profiles^{10,11,13,20} were used to deduce the composition and structure of the glycan species in each fraction. The results are summarized in Table 1 and Figure 2.

Optimization of the Labeling Step. The conditions for the labeling of oligosaccharides with APTS by reductive amination were optimized to obtain a robust, high-throughput microscale procedure with conservation of sialic acid residues. First, the preferred concentration of APTS for glycan labeling was determined: plasma samples were labeled using initial concentrations of 5, 10, 20, and 40 mM APTS. The results are depicted in Figure 3 showing an increasing yield of APTS-

labeled N-glycans with increasing APTS concentration. Because APTS is a rather expensive label, we chose to use 20 mM APTS for future experiments even though 40 mM APTS appeared to give slightly better yields of labeled glycans.

Second, the type and amount of reducing agent were optimized. Besides the widely used NaBH₃CN, we have applied 2-picoline borane for the labeling of N-glycans with APTS. We recently introduced 2-picoline borane for fluorescent labeling of oligosaccharides with 2-AA as well as 2-AB.¹⁴ Interestingly, we could show for APTS-labeling that variation of the NaBH₃CN concentration resulted in biased labeling of certain glycans, as is illustrated by the changes in ratio Peak 2/Peak 32 (ratio of nonfucosylated bisialylated Glycan 2 (H5N4S2) over nonfucosylated nonsialylated Glycan 32 ((H5N4), see Table 1) in Figure 3. Selective labeling was, however, not observed using 1, 0.5, 0.2, 0.1, and 0.05 M 2-picoline borane as the reducing agent (Figure 3). As the highest yields of labeled glycans were observed with addition of 0.2 M of 2-picoline borane (Figure




Figure 3. Optimization of the labeling procedure for labeling of N-glycans with the fluorescent tag APTS. Both the absolute intensity (in RFU) of the highest peak (Peak 2) as well as the ratio of nonfucosylated bisialylated Glycan 2 over nonfucosylated nonsialylated Glycan 32 are depicted.

3), it was decided to use 0.2 M of 2-picoline borane for the labeling of N-glycans by APTS to minimize bias and discrimination in the labeling step.

As glycan labeling by reductive amination is acid-catalyzed, the type and concentration of acid used in the reaction may influence the yield of APTS labeled glycans. This has indeed been reported previously, where stronger acids were associated with better yields.²² It is also known, however, that sialic acids can be hydrolyzed from the intact glycans under acidic conditions.¹⁷ Therefore, five different acids were applied as catalyst for the labeling of N-glycans with APTS: acetic acid ($pK_a = 4.74$), citric acid ($pK_a = 3.13, 4.76, 6.40$), malonic acid ($pK_a = 2.83, 5.13$), trichloroacetic acid (TCA, $pK_a = 0.70$) and trifluoroacetic acid (TFA, $pK_a = 0.30$), all at an initial concentration of 1.2 M. It was observed that citric acid resulted in the highest yields. Because the use of citric acid also resulted in only very modest loss of sialic acids as monitored by the ratio between sialylated glycans and nonsialylated glycans (Peak 2/Peak 32), it was decided to use citric acid in the following experiments. Citric acid was then applied at four different concentrations (see Figure 3). It was observed that acid hydrolysis was more pronounced at higher acid concentrations. However, the yield of labeled glycans (intensity of Peak 2) also increases with increased acid concentration (Figure 3). A initial citric acid concentration of 3.6 M was chosen for further experiments as it results in high yields and sialic acid residues are largely retained.

Optimization of the Purification Step. Fluorescent dye as well as reducing agent was employed in large excess to obtain efficient labeling. After derivatization, APTS-labeled glycans have to be purified from the reaction mixture, which contains large amounts of label, reducing agent, acid, salts, SDS, NP-40

as well as protein. Previously, the use of size exclusion has been reported for purification of APTS labeled N-glycans.^{10,11,23} However, these sample preparation methods resulted in either rather tedious procedures or largely impure samples, still containing substantial amounts of the derivatization agent. Recently, we reported a HILIC-based high-throughput purification method for 2-AA-labeled N-glycans. As the retention of labeled glycans in HILIC is mainly based on the glycan part, it was expected that HILIC would be a good choice for the purification of APTS-labeled glycans as well. Four different HILIC phases (cellulose, Sepharose, diol-bonded silica beads, and Biogel P-10 (a polyacrylamide based phase)) were examined for the purification of APTS labeled glycans, and results were compared to those obtained by the previously published, effective but labor intensive size exclusion chromatography (SEC) procedure¹¹ (see Figure S1, Supporting Information). Clearly, the samples purified using cellulose or the polyacrylamide (Biogel P-10) phase showed electropherograms that were very similar to those obtained after SEC purification. Since cellulose is a polysaccharide phase and bleeding may occur, the polyacrylamide phase was preferred and subsequent optimization of the purification procedure was performed using this stationary phase.

As a next step, the amount of the Biogel P-10 stationary phase was optimized. For all the tested amounts of Biogel P-10 (1, 2, 3, 5, and 10 mg) the ratio sialylated/nonsialylated (ratio Peak 2/Peak 32) remained constant, and with larger amounts of stationary phase, higher peak intensities for all glycan-peaks were obtained (data not shown). An amount of 10 mg of Biogel P-10 was therefore used for further optimization. Moreover, the washing conditions were optimized. Washing using acetonitrile/water (80:20 v/v) was tested, however, large amounts

of APTS could still be observed in the sample (data not shown). To elute APTS under washing conditions, we aimed to disrupt possible interactions of the strong negative charges of APTS with the stationary phase and therefore it was decided to introduce triethylamine (TEA) to the washing solution. Concentrations of TEA (400 mM, 250 mM and 100 mM) were added to the acetonitrile/water washing solution (80:20 v/v). As all three concentrations resulted in a similar reduction of APTS in the eluates (data not shown), the lowest TEA concentration of 100 mM TEA was used for further experiments. We then tested the effect of lowering the pH of the TEA-containing washing solution by addition of acetic acid. The pH of the washing solution containing 100 mM TEA was adjusted to 10 and 8.5 using acetic acid. The most efficient removal of APTS was observed at pH 8.5 (data not shown), and these conditions were selected for further experiments. To examine the number of washing steps necessary to clean the derivatized glycans from excess APTS, washing was performed three, five and seven times using 200 μ L solvent. Increasing the number of washing steps from three to five resulted in better purified samples, while a further increase to seven washing steps did not show any further reduction of APTS (data not shown). To avoid excess salts in the eluate due to the TEA/acetic acid in the washing solvent, which would most probably interfere with the electrokinetic injection in the CGE, three additional washing steps of 200 μ L acetonitrile/water (80:20, v/v) were performed. Prior to elution of the APTS-labeled glycans, a 100 μ L volume of water was added to the Biogel P-10 beads. The change from acetonitrile/water (80:20, v/v) to water made the beads swell. Elution was then achieved using 2×200 μ L water. Eluates were combined for further analysis. In a third elution step hardly any APTS-labeled glycans could be observed.

Sample Application. Initially, injection solutions for CGE-LIF were prepared diluting the samples 1:10 in 30 μ L water prior to analysis by CGE-LIF. However, total signal intensities of the APTS-glycan peaks showed relative standard deviations up to 50% for replicate injections. A possible explanation for such large variation could be different salt concentrations in the sample, which influence the electrokinetic injection. Therefore, in a first attempt, 1, 4, and 10% of running buffer (void of sieving matrix) were added to the sample. Signal intensity was largely decreased already using 1% of running buffer, while the variation was still very high (data not shown). Therefore, two other options were examined: the addition of Hi-Di Formamide was tested as suggested by Applied Biosystems for sequencing analysis. Moreover, DMSO was chosen as it is a good solvent for glycans and is less volatile than aqueous solvents preventing the drying of samples during overnight runs. Interestingly, dilution of the samples in DMSO resulted in largely enhanced signal intensities, as well as much lower variation between replicates. Optimal sample application conditions were achieved by adding 2 μ L of aqueous sample to 60 μ L of DMSO.

Repeatability of the Method. A 96-well plate containing 48 aliquots of a single plasma sample from a healthy donor, subjected to N-glycan release, APTS labeling and HILIC purification is called “batch” in the following. Three batches were prepared at different days. For determination of the measurement repeatability per capillary batch 1 was analyzed by CGE-LIF on three consecutive days. To this end, samples after HILIC purification were stored in the freezer and thawed on two days for second and third CGE-LIF runs. Aliquots were taken freshly thereof for each of the consecutive analyses. Notably, per sample the same capillary was used for the analysis on the three




Figure 4. Repeatability of plasma N-glycome analysis by CGE-LIF. (A) Average relative peak heights of the 12 major peaks are blotted together with the intrabatch standard deviation. (B) Relative standard deviations are given for the repetition per capillary and the intrabatch as well as the interbatch repeatability test.

days. Data were aligned, and peak heights were determined for 12 major peaks (Figure 4A). RSDs were determined by comparing the signals obtained from CGE-LIF analyses on three consecutive days for a specific plate position (well), thereby excluding the effect of variability between capillaries. The RSDs of the specific well positions were averaged over all 48 samples providing the intrabatch repeatability per capillary and are plotted in Figure 4B. RSD values were approximately 2% for the 7 early eluting peaks and approximately 5% for the 5 late eluting peaks (Figure 4B), resulting in an average RSD of 3.2% for the 12 peaks analyzed.

To assess the overall intrabatch and interbatch repeatability, two batches containing 48 replicates of the same single plasma sample were analyzed on different days involving glycan release, APTS labeling, HILIC, and multiplexed CGE-LIF. For intrabatch repeatability RSDs were determined per peak of the 48 CGE-LIF electropherograms obtained from one run. This was done for the three batches, and the results were averaged. RSD values were approximately 4% for the 7 early eluting peaks and approximately 9% for the 5 late eluting peaks (Figure 4B), resulting in an average RSD of 6.2% for the 12 peaks analyzed. For interbatch repeatability RSDs were calculated per peak for 3×48 electropherograms obtained for the three different




Figure 5. Changes in plasma glycosylation with different stages during and after pregnancy. Integrals of 12 different N-glycan peaks are plotted relative to the bisected bisialylated Glycan 2. Integrals were determined for samples taken at six different time points: in the 1st trimester of pregnancy (1), in the 2nd trimester of pregnancy (2), in the 3rd trimester of pregnancy (3), 6 weeks postpartum (4), 3 months postpartum (5), and 6 months postpartum (6). As the last time point is regarded "normal", the corresponding values were set to 100%. Time series are plotted for Glycans 1 (A), 4 (B), 5 (C), 6 (D), 14 (E), 17 (F), 20 (G), 24 (H), 28 (I), 30 (J), 32 (K), and 34 (L).

batches. Again, early eluting peaks showed a lower interbatch RSD than late-eluting peaks. The average interbatch RSD of the 12 selected peaks was found to be 15.8%.

Registration of Glycosylation Changes with Pregnancy. In a next step we wanted to test whether the CGE-LIF method was able to show glycosylation changes associated with a biological process. Sera from a healthy woman with an uncomplicated pregnancy were taken in the first, second and third trimester of her pregnancy, as well as 6 weeks, 3 months and 6 months postpartum. The total N-glycome was analyzed for all six sera. Alterations of IgG glycosylation patterns are known to be associated with different stages of pregnancy.^{15,24} As expected we could observe altered glycosylation patterns in total serum N-glycans, as is illustrated in Figure 5. Integrals of the 13 most abundant Peaks (1, 2, 4, 5, 6, 14, 17, 20, 24, 28, 30, 32, 34) were calculated and normalized to the biantennary bisialylated Glycan 2. The glycosylation status at 6 months postpartum was regarded 'normal', and therefore, the relative abundance of each glycan over the whole observation period was normalized to 100% at time point 6. Clearly, the relative levels of different glycans changed during and shortly after

pregnancy. Interestingly, increases up to 170% (Glycan 4, third trimester) as well as decreases down to 62% (Glycan 30, third trimester) could be observed. Moreover, differences in expression time courses were observed between biantennary sialylated glycans with bisecting GlcNAc (Glycans 6 and 20) and these species lacking bisecting GlcNAc (Glycans 5 and 17): the nonbisected species were elevated during pregnancy, while the bisected species were not. While these data were only obtained for one pregnant woman and will therefore not be suitable to deduce general patterns of plasma glycosylation changes during pregnancy, they do clearly show that the presented high-throughput method for analysis of APTS-glycan by CGE-LIF is suitable for detecting physiological glycosylation changes.

Discussion

In the field of -omics analysis, glycomics is increasingly being used to generate valuable information on the glycosylation of, for example, serum IgGs or total plasma proteins. Several recent studies have shown the potential of this approach to provide novel insight in unraveling stages of healthy vs diseased (e.g.,

refs 25–27). However, clinical studies, and certainly those studies that are performed at the population level, are characterized by large sample numbers, necessitating the application of fast and robust analysis procedures. Over the years, several fast procedures using direct mass spectrometry have been developed, though separation of isobaric structures is impossible using this approach. We and others recently published on fast and high-throughput sample preparation procedures for the labeling of N-glycans with 2-AA¹⁶ or 2-AB,⁵ however, analysis of these samples, other than by direct mass spectrometry is still tedious, as multiplexed LC systems are not readily available.

An alternative analysis approach to MS and LC-based techniques has been published^{9,10} in which a commercially available DNA-sequencer, which is commonly used in genetics departments, was used for the analysis of APTS-labeled N-glycans. In a typical “genomics” approach, 96 samples can be routinely analyzed in parallel, which indicates the potential of this method for high-throughput glycan profiling. Unfortunately, for glycan analysis applications, so far no satisfying sample preparation has been described, hampering the application of this methodology for large-scale studies. Indeed, publications so far have mainly described the use of 1–8 channel systems.^{10,11}

We here describe such a high-throughput sample preparation procedure for the facile generation of APTS-labeled N-glycans from human plasma proteins. The procedure is performed at the 96-well plate level, and the whole procedure can be performed with a hands-on time of 2.5 h. The protocol is suitable for automation, and performing 2 plates in parallel reduces net labor allocation time by at least 30 min per plate. As two overnight incubation steps are necessary, one for the release of the N-glycans and one for the labeling of the N-glycans, the total procedure takes 2.5 days. The costs for the preparation of APTS-labeled glycans at the 96-well plate level are comparable to those for the 2-AA labeling procedure described previously,¹⁶ though APTS is much more expensive than 2-AA, the major cost factor in sample preparation is the PNGase F enzyme.

Glycan labeling was optimized to preserve sialic acid residues. While others remove sialic acid residues in order to reduce sample complexity,^{9,10} we chose to monitor sialylation, because plasma protein sialylation is known to be involved in determining the protein half-life with loss of sialic acids making galactose residues accessible targets for protein removal by the asialoglycoprotein receptor.²⁸ Moreover, hypersialylation of plasma proteins has been described in cancer,²⁹ and IgG sialylation has been described as an anti-inflammatory signal.^{30,31}

The use of 2-picoline borane in the present protocol has two major advantages compared to the conventional use of NaBH₃CN. First, application of 2-picoline borane does not result in the release of HCN gas, and is thus less harmful for researchers and environment, especially when using high-throughput formats. Second, we showed that labeling of N-glycans with APTS using NaBH₃CN results in biased labeling of specific glycans with a strong dependence on the concentration of reducing agent (Figure 2). Such biased labeling was not observed for 2-picoline borane. Thus, the 2-picoline borane method is less susceptible to changes in the glycosylation pattern due to differences in concentration of reagents and, therefore, more robust.

We chose to test various HILIC stationary phases for purification of APTS-labeled N-glycans. Next to the polysac-

charide-based stationary phases cellulose and Sepharose, we tested a diol-bonded stationary phase for HILIC-SPE. Diol-bonded stationary phases have previously been described for carbohydrate analysis by HILIC-HPLC–ESI–MS.³² Moreover, the polyacrylamide (PAA) stationary phase Biogel P-10 was tested. While this material is commonly used for SEC, we were able to apply it in HILIC-SPE mode using solvent mixtures with 80% acetonitrile for binding and washing, and water for elution of APTS-labeled glycans. The use of a PAA-containing stationary phase for the separation of aminopyridine-labeled oligosaccharides, has been described before in HILIC-HPLC–ESI–MS.³³ Among the tested HILIC stationary phases, Biogel P-10 was chosen as it resulted in a CGE-LIF total plasma N-glycome profile which was very similar to the profile obtained after APTS-glycan purification by SEC (Figure S1, Supporting Information). The APTS-glycan profile obtained after SEC purification¹¹ served as a reference and is considered to be unbiased as this purification method (ideally) does not involve adsorption. The amount of Biogel P-10 HILIC-SPE material was varied (1–10 mg stationary phase). As expected, the overall amount of captured glycans varied, but again, we did not observe significant differences in the profiles and thus the relative intensities, demonstrating the robustness of the method for relative quantification (data not shown).

Using Biogel P-10, glycan samples cannot be contaminated with polysaccharides by bleeding of the stationary phase, as PAA is the only stationary phase in our test, which is not polysaccharide based. The HILIC-SPE procedure was optimized with regard to washing conditions. A washing solution containing 100 mM TEA/acetic acid (pH 8.5) in 80% acetonitrile was found to be largely efficient in eluting free APTS. Obviously, with the positively charged TEA, ionic interaction between APTS and the stationary phase are suppressed by ion pairing. Additional washing steps with 80% acetonitrile were included in the protocol in order to get rid of salt prior to elution, thus avoiding high and differing salt content in the injection solution for CGE. APTS-labeled glycans were then eluted with water and could directly be subjected to CGE-LIF analysis.

Electrokinetic injection on the CGE-LIF system was found to be particularly efficient when a small volume of aqueous HILIC-SPE eluate was diluted with DMSO. DMSO is a good solvent for carbohydrates and, apparently, also for APTS-labeled glycans. Moreover, it shows a low volatility, and sample volumes will thus hardly change by evaporation, which is particularly relevant when long series of measurements are performed on the DNA analyzer. The high peak intensities we obtained using a high percentage of DMSO in the sample injection solution might be explained by two facts: First, DMSO is a nonprotic solvent and the conductivity of the sample injection solution with a high percentage of DMSO is largely decreased. In addition, anions generally exhibit an increased pK_a value in DMSO compared to water.³⁴ This further reduces the conductivity of the sample injection solution as weak acids like acetic acid present from the sample workup do not account for conductivity. Only the very strong acids, like the sulfonic acids of the APTS label keep their high dissociation constants and thus exhibit the high electrophoretic mobility necessary for electrokinetic injection. It is interesting to note that the degree of sample dilution did not change the relative peak intensities, which shows, that there is no interference between the glycan species in quantification. This is in contrast to mass spectrometry, where quenching effects have to be taken into account.

The repeatability test showed that intrabatch RSDs for peak intensities are relatively low (6.2% on average) while the interbatch RSDs are considerably larger (15.8%). Higher RSDs are mainly observed for the late-eluting peaks which contain noncharged glycans. When large cohorts are being analyzed it may be necessary to address the interbatch variability by performing batch-corrections per plate for statistical analysis. This need is even stronger in our system, as the multiplexing of the analyses induces additional variation compared to single-capillary systems.

A major effort was made to annotate several of the N-glycan peaks in the electropherograms, as it was impossible to couple the separation system directly with mass spectrometry. Next to the comparison of migration positions with those of standard glycans, we relied on a HILIC-HPLC fractionation of a human plasma N-glycan sample. Aliquots of each fraction were labeled with 2-AA for high-resolution MALDI analysis and with APTS for CGE-LIF. Using this combined approach, a majority of 35 peaks was annotated. For future work, we will build up a database for the CGE-LIF electropherograms of APTS labeled N-glycans similar to Glycobase, a database for retention of 2-AB-labeled N-glycans on HILIC-HPLC columns.³⁵ The use of exoglycosidases, as well as more extensive, additional fractionation on a complementary separation system will however be necessary to improve the annotation. Glycans of identical composition were observed within distinct peaks (e.g., Peak 17 and Peak 21, see Table 1), demonstrating the separation of structural isomers. A similar separation was observed for the structural isomers Glycan 28 and Glycan 30, in accordance with literature.

The method was successfully applied to demonstrate changes in the total plasma N-glycome during pregnancy. Changes in the total plasma N-glycome may be caused by changes in the plasma concentration of specific glycoproteins. Alternatively, the glycosylation of an already present glycoprotein present at a high concentration may be changed. As total plasma was used, the identity of this possible protein cannot be addressed. An example for the latter is human plasma IgG, for which changes in Fc glycoforms have been described with age, pregnancy, and disease.^{15,24,36}

Conclusion

Labeling, purification and sample injection in CGE-LIF of APTS-derivatized human plasma N-glycans were optimized with respect to sensitivity, robustness, and minimized loss of sialic acids due to acid-hydrolysis. Peak identification of the majority of signals, also of overlapping signals, was achieved using HILIC-HPLC fractionation and MALDI analysis of glycans. We were able to show a high robustness of the method for its use in relative quantification. Peak identification of the majority of signals, including overlapping signals, was achieved using HILIC-HPLC fractionation and MALDI-MS analysis of glycans. The method was successfully applied to real samples showing biological variation.

Overall, we here present an N-glycan profiling method which can be performed with very short handling times. The method can be used for high-throughput glycosylation analysis. An advantage over "classical" HPLC or CE methods is that the analysis can be performed for up to 96 samples in parallel, with a very robust, easy to handle instrument, resulting in the separation of several structural isomers with high resolution.

Acknowledgment. We thank E. Steenvoorden for technical assistance. This work was supported by IOP grant IGE05007 and the Max Planck Society.

Supporting Information Available: Supplementary Figure showing the different electropherograms obtained using five different stationary phases. This material is available free of charge via the Internet at <http://pubs.acs.org>.

References

- (1) Apweiler, R.; Hermjakob, H.; Sharon, N. On the frequency of protein glycosylation, as deduced from analysis of the SWISS-PROT database. *Biochim. Biophys. Acta* **1999**, *1473* (1), 4–8.
- (2) Shental-Bechor, D.; Levy, Y. Folding of glycoproteins: toward understanding the biophysics of the glycosylation code. *Curr. Opin. Struct. Biol.* **2009**, *19* (5), 524–533.
- (3) Yu, Y. Q.; Gilar, M.; Kaska, J.; Gebler, J. C. A rapid sample preparation method for mass spectrometric characterization of N-linked glycans. *Rapid Commun. Mass Spectrom.* **2005**, *19* (16), 2331–2336.
- (4) Chu, C. S.; Ninonuevo, M. R.; Clowers, B. H.; Perkins, P. D.; An, H. J.; Yin, H.; Killeen, K.; Miyamoto, S.; Grimm, R.; Lebrilla, C. B. Profile of native N-linked glycan structures from human serum using high performance liquid chromatography on a microfluidic chip and time-of-flight mass spectrometry. *Proteomics* **2009**, *9* (7), 1939–1951.
- (5) Royle, L.; Campbell, M. P.; Radcliffe, C. M.; White, D. M.; Harvey, D. J.; Abrahams, J. L.; Kim, Y. G.; Henry, G. W.; Shadick, N. A.; Weinblatt, M. E.; Lee, D. M.; Rudd, P. M.; Dwek, R. A. HPLC-based analysis of serum N-glycans on a 96-well plate platform with dedicated database software. *Anal. Biochem.* **2008**, *376* (1), 1–12.
- (6) Ruhaak, L. R.; Zauner, G.; Huhn, C.; Bruggink, C.; Deelder, A. M.; Wuhler, M. Glycan labeling strategies and their use in identification and quantification. *Anal. Bioanal. Chem.* **2010**, *397* (8), 3457–3481.
- (7) Zhuang, Z.; Starkey, J. A.; Mechref, Y.; Novotny, M. V.; Jacobson, S. C. Electrophoretic analysis of N-glycans on microfluidic devices. *Anal. Chem.* **2007**, *79* (18), 7170–7175.
- (8) Bowman, M. J.; Zaia, J. Comparative glycomics using a tetraplex stable-isotope coded tag. *Anal. Chem.* **2010**, *82* (7), 3023–3031.
- (9) Callewaert, N.; Geysens, S.; Molemans, F.; Contreras, R. Ultrasensitive profiling and sequencing of N-linked oligosaccharides using standard DNA-sequencing equipment. *Glycobiology* **2001**, *11* (4), 275–281.
- (10) Laroy, W.; Contreras, R.; Callewaert, N. Glycome mapping on DNA sequencing equipment. *Nat. Protoc.* **2006**, *1* (1), 397–405.
- (11) Schwarzer, J.; Rapp, E.; Reichl, U. N-glycan analysis by CGE-LIF: profiling influenza A virus hemagglutinin N-glycosylation during vaccine production. *Electrophoresis* **2008**, *29* (20), 4203–4214.
- (12) Vanderschaeghe, D.; Festjens, N.; Delanghe, J.; Callewaert, N. Glycome profiling using modern glycomics technology: technical aspects and applications. *Biol. Chem.* **2010**, *391* (2–3), 149–161.
- (13) Vanderschaeghe, D.; Szekrenyes, A.; Wenz, C.; Gassmann, M.; Naik, N.; Bynum, M.; Yin, H.; Delanghe, J.; Guttman, A.; Callewaert, N. High-throughput profiling of the serum N-glycome on capillary electrophoresis microfluidics systems: toward clinical implementation of GlycoHepatoTest. *Anal. Chem.* **2010**, *82* (17), 7408–7415.
- (14) Ruhaak, L. R.; Steenvoorden, E.; Koelman, C. A. M.; Deelder, A. M.; Wuhler, M. 2-Picoline-borane: a non-toxic reducing agent for oligosaccharide labeling by reductive amination. *Proteomics* **2010**, *10* (12), 2330–2336.
- (15) van de Geijn, F. E.; Wuhler, M.; Selman, M. H.; Willemsen, S. P.; de Man, Y. A.; Deelder, A. M.; Hazes, J. M.; Dolhain, R. J. Immunoglobulin G galactosylation and sialylation are associated with pregnancy-induced improvement of rheumatoid arthritis and the postpartum flare: results from a large prospective cohort study. *Arthritis Res. Ther.* **2009**, *11* (6), R193.
- (16) Ruhaak, L. R.; Huhn, C.; Waterreus, W. J.; de Boer, A. R.; Neuss, C.; Hokke, C. H.; Deelder, A. M.; Wuhler, M. Hydrophilic interaction chromatography-based high-throughput sample preparation method for N-glycan analysis from total human plasma glycoproteins. *Anal. Chem.* **2008**, *80* (15), 6119–6126.
- (17) Bigge, J. C.; Patel, T. P.; Bruce, J. A.; Goulding, P. N.; Charles, S. M.; Parekh, R. B. Nonselective and efficient fluorescent labeling of glycans using 2-amino benzamide and anthranilic acid. *Anal. Biochem.* **1995**, *230* (2), 229–238.
- (18) Packer, N. H.; Lawson, M. A.; Jardine, D. R.; Redmond, J. W. A general approach to desalting oligosaccharides released from glycoproteins. *Glycoconj. J.* **1998**, *15* (8), 737–747.

- (19) Selman, M. H.; McDonnell, L. A.; Palmblad, M.; Ruhaak, L. R.; Deelder, A. M.; Wuhler, M. Immunoglobulin G glycopeptide profiling by matrix-assisted laser desorption ionization Fourier transform ion cyclotron resonance mass spectrometry. *Anal. Chem.* **2010**, *82* (3), 1073–1081.
- (20) Vanderschaeghe, D.; Debruyne, E.; Van Vlierberghe, H.; Callewaert, N.; Delanghe, J. Analysis of gamma-globulin mobility on routine clinical CE equipment: Exploring its molecular basis and potential clinical utility. *Electrophoresis* **2009**, *30* (15), 2617–2623.
- (21) Stumpo, K. A.; Reinhold, V. N. The N-Glycome of Human Plasma. *J. Proteome Res.* **2010**, *9* (9), 4823–4830.
- (22) Evangelista, R. A.; Guttman, A.; Chen, F. T. Acid-catalyzed reductive amination of aldoses with 8-aminopyrene-1,3,6-trisulfonate. *Electrophoresis* **1996**, *17* (2), 347–351.
- (23) Vanderschaeghe, D.; Laroy, W.; Sablon, E.; Halfon, P.; Van Hecke, A.; Delanghe, J.; Callewaert, N. GlycoFibroTest: DNA sequencer-based serum protein glycomics yields a highly performant liver fibrosis biomarker. *Mol. Cell. Proteomics* **2009**, *5*, 986–949.
- (24) Rook, G. A. W.; Steele, J. F.; Brealey, R.; Whyte, A.; Isenberg, D.; Sumar, N.; Nelson, J. L.; Bodman, K. B.; Young, A.; Roitt, I. M.; Williams, P.; Scragg, I. G.; Edge, C. I.; Arkwright, P. D.; Ashford, D.; Wormald, M. R.; Rudd, P. M.; Redman, C. W. G.; Dwek, R. A.; Rademacher, T. W. Changes in IgG glycoform levels are associated with remission of arthritis during pregnancy. *J. Autoimmun.* **1991**, *4*, 779–794.
- (25) Saldo, R.; Royle, L.; Radcliffe, C. M.; Abd Hamid, U. M.; Evans, R.; Arnold, J. N.; Banks, R. E.; Hutson, R.; Harvey, D. J.; Antrobus, R.; Petrescu, S. M.; Dwek, R. A.; Rudd, P. M. Ovarian cancer is associated with changes in glycosylation in both acute-phase proteins and IgG. *Glycobiology* **2007**, *17* (12), 1344–1356.
- (26) Ruhaak, L. R.; Uh, H. W.; Beekman, M.; Koeleman, C. A. M.; Hokke, C. H.; Westendorp, R. G.; Wuhler, M.; Houwing-Duistermaat, J. J.; Slagboom, P. E.; Deelder, A. M. Decreased levels of bisecting GlcNAc glycoforms of IgG are associated with human longevity. *PLoS One* **2010**, *5* (9), e12566.
- (27) Callewaert, N.; Van, V. H.; Van, H. A.; Laroy, W.; Delanghe, J.; Contreras, R. Noninvasive diagnosis of liver cirrhosis using DNA sequencer-based total serum protein glycomics. *Nat. Med.* **2004**, *10* (4), 429–434.
- (28) Stockert, R. J. The asialoglycoprotein receptor: relationships between structure, function, and expression. *Physiol. Rev.* **1995**, *75* (3), 591–609.
- (29) Arnold, J. N.; Saldo, R.; Hamid, U. M.; Rudd, P. M. Evaluation of the serum N-linked glycome for the diagnosis of cancer and chronic inflammation. *Proteomics* **2008**, *8* (16), 3284–3293.
- (30) Anthony, R. M.; Wermeling, F.; Karlsson, M. C.; Ravetch, J. V. Identification of a receptor required for the anti-inflammatory activity of IVIG. *Proc. Natl. Acad. Sci. U.S.A.* **2008**, *105* (50), 19571–19578.
- (31) Kaneko, Y.; Nimmerjahn, F.; Ravetch, J. V. Anti-inflammatory activity of immunoglobulin G resulting from Fc sialylation. *Science* **2006**, *313* (5787), 670–673.
- (32) Tolstikov, V. V.; Fiehn, O. Analysis of highly polar compounds of plant origin: combination of hydrophilic interaction chromatography and electrospray ion trap mass spectrometry. *Anal. Biochem.* **2002**, *301* (2), 298–307.
- (33) Ikegami, T.; Horie, K.; Saad, N.; Hosoya, K.; Fiehn, O.; Tanaka, N. Highly efficient analysis of underivatized carbohydrates using monolithic-silica-based capillary hydrophilic interaction (HILIC) HPLC. *Anal. Bioanal. Chem.* **2008**, *391* (7), 2533–2542.
- (34) Sarmini, K.; Kennedler, E. Ionization constants of weak acids and bases in organic solvents. *J. Biochem. Biophys. Methods* **1999**, *38* (2), 123–137.
- (35) Campbell, M. P.; Royle, L.; Radcliffe, C. M.; Dwek, R. A.; Rudd, P. M. GlycoBase and autoGU: tools for HPLC-based glycan analysis. *Bioinformatics* **2008**, *24* (9), 1214–1216.
- (36) Arnold, J. N.; Wormald, M. R.; Sim, R. B.; Rudd, P. M.; Dwek, R. A. The impact of glycosylation on the biological function and structure of human immunoglobulins. *Annu. Rev. Immunol.* **2007**, *25*, 21–50.

PR100802F